


GLUCOMMANDER DOWNTIME

If the Glucommander system has a downtime you can print the Glucommander downtime form from two places:

- 1. Desktop ICON (G)
- KnowledgLink → Policies and Guidelines → Downtime
 Procedures → Glucommander for IV Insulin





Recommended Procedure

- Obtain blood glucose 10 minutes prior to downtime
- Adjust insulin infusion per Glucommander recommendations
- Document current multiplier
- If Glucommander downtime is less than 1 hour resume using next blood glucose check
- If Glucommander downtime is greater than 1 hour use the Glucommander downtime form accessed per above to determine hourly titrations
- When Glucommander is back up, access Glucommander via the SmartClick link on the Flowsheet or MAR and restart Glucommander with the last multiplier from the downtime form (see tip sheet for editing IV multiplier after downtime)



