

Protocol Definition

Dräger RS 232 MEDIBUS

WARNING!

Strictly follow the Instructions for Use.
Any use of the software protocol
requires full understanding and strict
observation of these instructions. The
software protocol is only to be used for
purposes specified here.

Revision Level 6.00

Contents

Intended Use	3	Configure Realtime Transmission Command .	34
		Realtime Configuration Changed Command .	34
MEDIBUS Introduction	4	MEDIBUS Realtime-Extension Responses .	35
Initializing Communication	4	Realtime Configuration Response	35
Terminating Communication	5	MEDIBUS Realtime-Extension Realtime-Data	
Time-Out	5	Records	36
Allowable Characters	5	Structure of Realtime-Data Records	36
Software Handshaking	6	Sync-Commands	38
Suspending data transmission	6	Appendix	40
Resuming data transmission	6	MEDIBUS-Life-Cycle	40
Aborting data transmission	6	MEDIBUS Realtime-Extension Life-Cycle	42
MEDIBUS Commands	7	ASCII HEX Format	45
Structure of Commands	7	MEDIBUS example 1	46
Command Codes	8	MEDIBUS example 2	47
Request Trend Data Command	11	MEDIBUS example 3	48
MEDIBUS Responses	12	MEDIBUS example 4	49
Structure of Responses	12	MEDIBUS example 5	50
Responding to Control Commands and unknown		MEDIBUS example 6	51
Commands	13	Frequently Asked Questions	55
Responding to corrupt Commands	13	Logbook of Changes	58
Responding to Data Request Commands	14	Table of used ASCII-Codes	59
MEDIBUS Device Extension for			
Infusion Pumps	22		
Introduction	22		
MEDIBUS Device Extension Protocol Concept	22		
Infusion Pump Commands	23		
Infusion Pump Communication Life Cycle	24		
Request Infusion Pump Configuration	25		
Infusion Pump Configuration Changed			
Command	26		
Infusion Pump Data Update Commands	27		
MEDIBUS Realtime-Extension Introduction	30		
Time-Out	30		
Allowable Characters for Realtime			
Transmission	31		
MEDIBUS Realtime-Extension Commands .	33		
Command Codes	33		
Request Realtime Configuration Command . . .	33		

Intended Use

MEDIBUS is a software protocol intended to be used for exchanging data between a Dräger medical device and external medical or non-medical devices via RS 232 interfaces.

WARNING

Data transferred via MEDIBUS interfaces are for information only and are not intended as a basis for diagnosis or therapy decisions.

MEDIBUS consists of two independent software protocols one for the transmission of "slow" and one for the transmission of "fast" data.

"Slow" data:

Generated or updated in intervals of the magnitude of seconds. This part is called MEDIBUS.

"Fast" data:

Intended for the transmission of e. g. realtime curves. This part of MEDIBUS is called the "Real-time-Extension".

This manual contains a general description of the protocol including formats of commands and responses.

For device dependent descriptions of supported commands and data sets, port hardware and configurations for Dräger devices please refer to the individual device-specific MEDIBUS documentation.

MEDIBUS Introduction

The MEDIBUS protocol distinguishes two basic types of messages:

- commands
- responses.

A command is transmitted by one device to request data from the other device or to control its function.

A response is transmitted by one device upon receipt of a command from the other device.

Responses may contain embedded commands.

Initializing Communication

(refer to MEDIBUS-life-cycle-diagram, see page 40)

- To initialize communication or to restart communication after a time-out, a device must send the "Initialize Communications Command" **ICC**. Refer to section "Control commands" for the format of commands.
- A device considers communication initialized after having received either a response to a transmitted ICC or an ICC from the other device. Refer to section "Responding to commands" for the format of responses.
- Commands embedded in a response to an ICC are disregarded.

Terminating Communication

- To stop the communication the "STOP" command has to be sent.
- The command echo has to be checked to make sure the "Stop Communication" command has been received correctly. Further commands from the linked device may be ignored until the communication has been reinitialized.

Time-Out

- Any pause in the data flow exceeding 3 seconds leads to a time-out, terminating the communication link. To resume communication after a time-out, the device must re-initialize communication by transmitting an ICC command (see section "Initializing communication", page 4). Whenever a device receives an ICC command, it must send a response to it (see section "Responses", page 12).
- After receiving a command the device has to send a complete response within 10 seconds.
- If there is no need for sending commands or responses the "NOP"-command has to be sent in 2-second-intervals to keep the communication alive.

Allowable Characters


- Printable ASCII characters.
- Control characters defined in this Instructions for Use.

Software Handshaking

Some control characters can be sent at any time to control the flow of data. They do not require responses.


Suspending data transmission

- If a receiving device wants the transmitting device to suspend transmission, the ASCII "DC1" character (11H) must be sent.
- Upon receipt of this character, the sending device will suspend any transmission immediately until it receives the ASCII "DC3" character.




Resuming data transmission

- To request the transmitting device to resume data transmission, the ASCII "DC3" character (13H) must be sent within 3 seconds. Else communication will be reinialized with an ICCcommand.



Aborting data transmission

- To request the other device to abort sending a response or a command, the ASCII "CAN" character (18H) must be sent.
- Upon receipt of this character, any transmission in progress will be immediately aborted. Communication may be restarted by repetition of the last sent command immediately.





MEDIBUS Commands

Structure of Commands

- A command is a string of ASCII characters transmitted by one device to request data from the other device or to control its function.
- A command may be embedded in the response to another command, but a new command must not be transmitted until the response to the previous command has been received.
- If, however, the response to a command has not been received in full within 10 seconds since the transmission of the last command byte, the command may be repeated or a new command may be transmitted.

Commands have one of the following formats:



ESC ASCII "escape" character (1BH)

Command-Code Single byte code specifying the command.

ARGUMENT The argument string is of variable length n, but n must not exceed 251(0FBH) bytes. The string consists of printable ASCII characters, either text or ASCII HEX numbers. The format of the different arguments are specified in the following sections.

CHECKSUM Least significant 8-bit sum of all preceding bytes beginning with "ESC" in ASCII HEX format (see section "ASCII HEX Format", page 45).

CR ASCII "carriage return" character (0DH)

NOTE

To avoid communication breakdown, commands with arguments must not be sent to devices with a MEDIBUS version less than 3.00.

Command Codes

Control Commands

Control commands are used to initialize, control and stop communication.

Command	Code
No Operation (NOP)	30H
Initialize Communication (ICC)	51H
Stop Communication (STOP)	55H

Data Request Commands

Data request commands are used to request data.

Command	Code
Request current Alarms (Codepage 3)	23H
Request current measured Data (Codepage 1)	24H
Request current low Alarm Limits (Codepage 1)	25H
Request current high Alarm Limits (Codepage 1)	26H
Request current Alarms (Codepage 1)	27H
Request current Date and Time	28H
Request current Device Setting	29H
Request current Text Messages	2AH
Request current measured Data (Codepage 2)	2BH
Request current low Alarm Limits (Codepage 2)	2CH
Request current high Alarm Limits (Codepage 2)	2DH
Request current Alarms (Codepage 2)	2EH
Request Device Identification	52H
Request Trend Data Status	6CH
Request Trend Data	6DH

Miscellaneous Commands

Command	Code
Time changed	49H
Configure Data Response Command	4AH

Time Changed Command

The "Time changed" command is sent if during runtime the time or date of the device has been changed. The receiving device can now ask for current date and time. Because of the "Time changed" command there is no need for sending periodically the "Request Current Date and Time" command.

Configure Data Response Command

The "Configure Data Response" command is used to limit the number of data responded by a data source on a data request command. On receipt of any of these data request commands the receiver has to send all actual valid data. In cases where only a few of the possible data are used, the requesting device may configure the responding device to send only these used data by sending a "Configure Data Response" command. The codes of useful data are given in the argument as follows:

DATA TYPE	DATA CODE 1	DATA CODE 2	...	DATA CODE n	
0	+1	+3	+5	+2n-2	+2n n <= 125

DATA TYPE:

One byte identifying the data type to configure.

This may be:

- 24H for current Data, low Alarm Limits and high Alarm Limits (codepage 1)
- 27H for current Alarms (codepage 1)
- 29H for current Device Settings
- 2AH for current Textmessages
- 2BH for current Data, low Alarm Limits and high Alarm Limits (codepage 2)
- 2EH for current Alarms (codepage 2)

DATA CODE:

Two byte ASCII HEX data code. See appendices for code numbers.

The configuration stays valid until receipt of a new "Configure Data Response" command. After re-initialization of communication (ICC) and if the "Configure Data Response" is send without data codes, the configuration is set to its default state, where internal programmed configuration is used.

An example is given in Appendix.

Request Trend Data Command

The "Request Trend Data" command (6DH) is used to request a sequence of samples of one trend parameter.

The command argument has the following format:

CODE PAGE	DATA CODE	COUNT	BEGIN	
0	+2	+4	+6	+14

CODE PAGE: Two byte ASCII HEX number identifying the code page of the trend parameter:

24H for data code page 1

2BH for data code page 2

DATA CODE: Two byte ASCII HEX number containing the data code of the trend parameter. A parameter's data code is the same as used for the "Request current measured data" response.

COUNT: Two byte ASCII HEX number defining the maximal number of trend samples allowed in the response.

BEGIN: Eight byte ASCII HEX number specifying the time stamp of the eldest requested trend sample. (Note: BEGIN is undefined if COUNT is 0.)

The time stamp is packed in accordance to the Microsoft Windows™ 32-bit date and time format:

NOTE

BEGIN is undefined if COUNT is 0.

YYYYYYMMMMDDDDhhhhmmmmsssss
MSB LSB

with

YYYYYY year relative to 1980 (range 0 - 119)

MMMM month (range 1 - 12)

DDDD day (range 1 - 31)

hhhh hour (range 0 - 23)

mmmmmm minute (range 0 - 59)

sssss second (range 0 - 29 in 2-second intervals)


NOTE

Prior to the "Trend Data Request" the device should request the trend data status in order to know the availability of trend data (refer to chapter "Trend Data Status Response" for the respective response).

MEDIBUS Responses

Structure of Responses

Upon receipt of a command, a device must respond to it within 10 seconds. A command may be embedded within the response. The following format has to be used:



The response is of variable length, but must not exceed 3845 bytes.

SOH	ASCII "Start of Header" character (01H)
Command ECHO	Echo of the command code being responded to.
RESPONSE	Data as requested by the command, see sections "MEDIBUS Specification"
CHECKSUM	Least significant 8-bit sum of all preceding bytes beginning with "SOH" in ASCII HEX format (see section "ASCII HEX Format", page 45).
CR	ASCII "carriage return" character (0DH)

Responding to Control Commands and unknown Commands

A response to a control command or unknown command acknowledges receipt of the command, but contains no data:

SOH	Command-ECHO	CHECKSUM	CR
0	1	2	4
5 Byte			

SOH ASCII "Start of Header" character (01H).

Command ECHO Echo of control command being responded to.

CHECKSUM Least significant 8-bit sum of all preceding bytes beginning with "SOH" in ASCII HEX format (see section "ASCII Hex format", see page 45).

CR ASCII "carriage return" character (0DH).

- Refer to section "Control commands" for currently defined control commands.

Responding to corrupt Commands

If the received command is corrupt (bad checksum), the command echo field must consist of an ASCII "NAK" character (15H) and there must be no response field:


SOH	NAK	CHECKSUM	CR
01H	15H	31H 36H	0DH
0	1	2	3
4			
5 Byte			

Responding to Data Request Commands

Responses to data request commands contain the current values of a device. In case a value isn't available at a certain time (for example caused by temporary measurement problems) nothing must be sent for this value. Vice versa this means, that a value is invalid if it is not included in a data request response

Current Measured Data and Alarm Limit Response

This response must be sent in reply to the "Request current measured Data (codepage 1)" command (24H), "Request current low Alarm Limits (codepage 1)" command (25H), "Request current high Alarm Limits (codepage 1)" command (26H), "Request current measured Data (codepage 2)" command (2BH), "Request current low Alarm Limits (codepage 2)" command (2CH) or "Request current high Alarm Limits (codepage 2)" command (2DH). It contains the current values of all measured parameters or alarm limits available on the responding device. The response field has the following format:



DATA CODE

Two byte ASCII HEX number identifying the parameter or alarm limit.

DATA

Four byte ASCII field containing the current value of the parameter or alarm limit. See appendix for data formats. Surplus character positions and leading zeros must be replaced by an ASCII "Space" (20H).

Alarm Status Response

This response must be sent in reply to the "Request current Alarms (codepage 1)" command (27H) or "Request current Alarms (codepage 2)" command (2EH). It contains the alarm priority, alarm code, and alarm message for all currently active alarms on the responding device. The response field has the following format:

ALARM Priority 1	ALARM CODE 1	ALARM PHRASE 1	ALARM Priority 2	ALARM CODE 2	ALARM PHRASE 2	---	ALARM Priority n	ALARM CODE n	ALARM PHRASE n	
0	1	3	15	16	18	30	n-15	n-14	n-12	n

ALARM Priority

One byte field specifying the alarm priority (number in the range of 1 to 31), 31 being the highest priority. The priority is encoded by adding 30H. The priorities, therefore, lie in the range from ASCII "1" (31H) to ASCII character "O" (4FH).

ALARM CODE

Two byte ASCII HEX number identifying the alarm.

ALARM PHRASE

Twelve byte ASCII character string describing the alarm.

Time & Date Update Response

This response is sent in reply to the "Request Current Date & Time" command (28H). It contains the current date and time from the responding device. The response field has the following format:

TIME (HH:MM:SS)	DATE (DD-MMM-YY)
0	+8
	+17
	017

TIME: Eight byte field containing ASCII numeric characters representing the current time in hours (HH), minutes (MM), and seconds (SS). Leading zeroes shall not be suppressed.¹⁾

DATE: Nine byte field containing ASCII alpha-numeric characters representing the current day (DD), month (MMM), and year (YY). The first three letters for each month are sent in ASCII. Leading zeroes shall not be suppressed.¹⁾

1) The PM 8040 substitutes leading zeros in the hours, day and year field with ASCII spaces (20H), e.g. ' 8:06:05 4-MAR 2'.

NOTE


The month representing letters (MMM) has to be sent in German language.

German month representation in time & Date update responses

MONTH	REPRESENTATIVE
January	JAN
February	FEB
March	MAR
April	APR
May	MAI
June	JUN
July	JUL
August	AUG
September	SEP
October	OKT
November	NOV
December	DEZ

Device Setting Responses

This response must be sent in reply to the "Request Current Device Settings" command (29H). It contains the current values of all device settings applicable with the responding device. The response field has the following format:




SETTING CODE Two Byte ASCII HEX number identifying the parameter.

SETTING Five byte ASCII field containing the current value of the specified parameter. Refer to appendix for the specific formats of the parameters. Surplus character positions and leading zeros must be filled up with ASCII "SPACE" (20H).

Text Message Response

This response must be sent in reply to the "Request Text Messages" command (2AH). It contains all the text messages the requested device currently holds for user information, along with the text code, text length and an end-of-text marker. The response field has the following format:



- TEXT CODE Two byte ASCII HEX number identifying the text messages.
- LENGTH One byte field specifying the text length, a number in the range from 1 to 32. A text must not be longer than 32 characters. The length is encoded to ASCII format by adding 30H to the decimal length value. Thus, the text length ranges from ASCII "1" (31H) to ASCII "P" (50H).
- TEXT ASCII character string. Refer to appendix for text messages.
- ETX End-of-text marker (ASCII-Code 03H).

NOTE
 The length of the response field is limited to 3840 bytes. Due to different lengths of text messages, the maximum number of text messages in a response field depends on the length of the individual text messages.

Device Identification Response

This response must be sent in reply to the "Request Device Identification" command (52H). It contains the identification number, name and release number of the responding device and the MEDIBUS release number. The response field has the following format:

ID NUMBER NNNN	NAME 'Device Name'	REVISION DD.DD:MM.MM
0	4	n+6
n+17 Byte		
0 < n ≤ 32		

020

ID NUMBER Four byte field containing the ASCII device identification number NNNN.

NAME ASCII character string delimited by apostrophes (ASCII Code 27H). Therefore, the device name itself must not contain apostrophes. The length of the device name may range from 1 to 32 characters.

REVISION Eleven byte field containing ASCII characters representing the device revision level (DD.DD) and the MEDIBUS revision level (MM.MM).

**All identification numbers will be defined by
Dräger Medical AG & Co. KG, Lübeck.**

Trend Data Status Response

This response sends information about the available trend samples for each trend parameter. It is sent on reply to the "Request Trend Data Status" command (6CH). The response field has the following format:

PARAMETER 1				PARAMETER 2				
CODE PAGE	DATA CODE	COUNT	BEGIN	CODE PAGE	DATA CODE	COUNT	BEGIN	
0	+2	+4	+10	+18	+20	+22	+28	+36

CODE PAGE: Two byte ASCII HEX number identifying the code page of the trend parameter:
 24H for data code page 1
 2BH for data code page 2

DATA CODE: Two byte ASCII HEX number containing the data code of the trend parameter. A parameter's data code is the same as used for the "Request current measured data" response.

COUNT: Six byte ASCII HEX number defining the number of trend samples available for the respective parameter.

Leading zeros may be replaced by ASCII spaces (20H).

BEGIN: Eight byte ASCII HEX number containing the time stamp of the eldest available trend sample of the respective parameter.

For the format of the time stamp refer to chapter "Request Trend Data Command".

Trend Data Response

This response is sent in reply to the "Request Trend Data" command (6DH). It contains a sequence of samples of one parameter with the respective time stamps.

The response field has the following format:

CODE PAGE	DATA CODE	COUNT = N	VALUE 1	TIME 1	VALUE 2	TIME 2	...	VALUE N	TIME N		N ≤ 255
0	+2	+4	+6	+10	+18	+22	+30	+n-12	+n-8	+n	n < 3066

- CODE PAGE:** Two byte ASCII HEX number identifying the code page of the trend parameter:
 24H for data code page 1
 2BH for data code page 2
- DATA CODE:** Two byte ASCII HEX number containing the data code of the trend parameter. A parameter's data code is the same as used for the "Request current measured data" response.
- COUNT:** Two byte ASCII HEX number defining the number of trend samples contained in the response.
- VALUE:** Four byte ASCII field containing the value of the trend sample. The format for each parameter is the same as used for the "Request current measured data" response.
- TIME:** Eight byte ASCII HEX number containing the time stamp of the trend sample.
 For the format of the time stamp refer to chapter "Request Trend Data Command".

MEDIBUS Device Extension for Infusion Pumps

Introduction

The standard MEDIBUS protocol is designed for data exchange between two medical devices, where each of the devices can support only one source of a certain data item. In opposite to the standard protocol the MEDIBUS device extension for infusion pumps is meant to support several pumps on one communication port, all pumps sourcing identical data items, e.g. the same measured data, alarms, settings etc.

MEDIBUS Device Extension Protocol Concept

The basic concept of the MEDIBUS device extension is stacking the MEDIBUS Device Extension command 59H with sub-command codes and addressing a specific pump by a channel number.

Stacking of command codes means that a device extension telegram always uses the leading command code 59H to indicate that the telegram is a MEDIBUS Device Extension command or response for infusion pumps. The leading command is followed by a sub-command code, e.g. specifying the actual data update request. If applicable the sub-command is followed by a channel number, defining to which device the telegram is related.


Infusion Pump Commands

Leading Command	Code
Device Extension Command for Infusion Pumps	59H

Sub-Commands	Code
Request current measured data	24H
Request current low alarm limits	25H
Request current high alarm limits	26H
Request current alarms	27H
Request current device settings	29H
Request current text messages	2AH
Request infusion pump configuration	41H
Infusion pump configuration changed	42H

Infusion Pump Communication Life Cycle

The following diagram shows a typical communication life cycle for a device being connected to an infusion pump device.



030

Request Infusion Pump Configuration

After communication has initialized the requesting device first has to ask the pump device for available pumps using the Request Infusion Pump Configuration sub-command (41H).

The pump device will respond with a list of identifier for available pump channels.

Later these pump channel identifier have to be used by the requesting device to request data from a specific pump.

Request Infusion Pump Configuration Command:

ESC	Lead-Command	Sub-Command 41H		Checksum D9H		CR
1BH	59H	34H	31H	44H	39H	0DH

Infusion Pump Configuration Response:

ESC	Lead-Command	Sub-Command 41H		Pump channel 01H		Pump channel xyH		Checksum		CR
01H	59H	34H	31H	30H	31H	xxH	yyH	CS1	CS2	0DH

Sub-Command: 2 bytes ASCII-Hex field.

Pump channel: 2 bytes ASCII-Hex field specifying the identifier of available infusion pump channels within the range of 0 to 255 (00H to FFH).

Infusion Pump Configuration Changed Command

If an infusion pump has been added to or removed from the pump device and the pump device has been asked for its configuration before, the pump device will issue an Infusion Pump Configuration Changed sub-command (42H).

After receiving the Infusion Pump Configuration Changed sub-command the requesting device has to restart the infusion pump communication again by sending the Request Infusion Pump Configuration sub-command.

Infusion Pump Configuration Changed Command:

ESC	Lead-Command	Sub-Command 42H		Checksum DAH		CR
1BH	59H	34H	32H	44H	41H	0DH

Infusion Pump Configuration Changed Response:

SOH	Lead-Command	Sub-Command 42H		Checksum C0H		CR
01H	59H	34H	32H	43H	30H	0DH

Sub-Command: 2 bytes ASCII-Hex field.

Infusion Pump Data Update Commands

Infusion pump data update commands are sent to request all actual values of a certain data type from a specific pump.

Infusion Pump Data Update Request:

To request actual measured data, alarm limits, alarms, settings or text messages the requesting device has to send the according data update sub-command (24H, 25H, 26H, 27H, 29H, 2AH - refer to chapter Infusion Pump Commands for the respective meaning of the sub-commands) and the according infusion pump channel identifier.

ESC	Lead-Command	Data Update Request Sub-Command	Pump Channel xyH	Checksum		CR
1BH	59H	COMMAND2	CHANNEL	CS1	CS1	0DH

COMMAND2: 2 bytes ASCII-Hex field specifying the requested data type (24H, 25H, 26H, 27H, 29H or 2AH).

CHANNEL: 2 bytes ASCII-Hex field specifying the identifier of the infusion pump channel.

Infusion Pump Data Update Response:

Responses to Infusion Pump Data Update commands contain all current values of one pump as specified by its channel identifier. In case a value isn't available at a certain time (for example caused by temporary measurement problems) nothing must be sent for this value. Vice versa this means, that a value is invalid if it is not included in a data update response

SOH	Lead-Command	Data Update Request Sub-Command	Pump Channel xy	RESPONSE FIELD	Checksum		CR
01H	59H	COMMAND2	CHANNEL	DATA 1 ... DATA n	CS1	CS1	0DH

COMMAND2: 2 bytes ASCII-Hex field specifying the responded data type (24H, 25H, 26H, 27H, 29H or 2AH).

CHANNEL: 2 bytes ASCII-Hex field specifying the identifier of the infusion pump channel.

DATA n: For alarms and text messages the response field contains a list of data with the same format as defined before for the according data type in the standard MEDIBUS protocol.

NOTE

For measured data, low and high alarm limits and settings the format of the response field deviates from the according format in the standard MEDIBUS protocol. Instead of 4 or 5 bytes in the standard protocol, for infusion pumps these values have a length of 8 bytes.

Example:

To request actual settings from pump channel 2:

ESC	Lead-Command	Setting Update Request Sub-Command 29H		Pump Channel 02H		Checksum 41H		CR
1BH	59H	32H	39H	30H	32H	34H	31H	0DH

The setting response for pump on channel 2 with two settings at code 01H with value '1234.5' and at code 03H with value '5' would be as follows:

SOH	Lead-Command	Data Update Request Sub-Command 29H		Pump Channel 02H		RESPONSE FIELD		Checksum 6DH		CR
01H	59H	32H	39H	30H	32H	DATA 1	DATA 2	36H	44H	0DH

With the RESPONSE FIELD as follows:

DATA 1										DATA 2									
Code 1		Value 1								Code 2		Value 2							
01H		1234.5								03H		5							
30H	31H	20H	20H	31H	32H	33H	34H	2EH	35H	30H	33H	20H	20H	20H	20H	20H	20H	20H	35H

MEDIBUS Realtime-Extension Introduction

The Realtime-Extension of MEDIBUS is designed to run independently and without any interference with the slow MEDIBUS concurrently on the same port. Since realtime data are usually created and transmitted at fixed time intervals, it may happen that slow and fast records are mixed. This means that transmission of a slow command or response may be interrupted by a record of fast data at any time.

Time-Out

Realtime bytes (Sync-Byte, Sync-commands, Realtime values) do not affect the 3 sec time-out described in the "Time-out" chapter of the "slow" MEDIBUS-Protocol.


Allowable Characters for Realtime Transmission

Slow communication uses ASCII codes up to 7FH, meaning the most significant bit (bit 7) is always 0. Realtime ("fast") data is distinguished from slow data in that the most significant bit (realtime data flag) is set. All received data < 80H must be ignored within realtime transmission.

There are four types of data-bytes which can be distinguished unequivocal:


Command Response Bytes

Slow communication's Control-Commands and the responses they require are transmitted in the following form:




Sync Byte

A Sync Byte is always the first byte in a realtime-databyte (see page 36). A Sync Byte is transmitted in the following format:




Sync-Command Bytes

Sync-Commands are used for information which have to be transmitted without delay (see section Sync-Commands, page 38). The form of a sync-command-byte and its argument is as follows:



Realtime-Value-Bytes

A realtime value is transmitted at a resolution of 12 bits. For transmission the value is divided into two data bytes:



MEDIBUS Realtime-Extension Commands

In addition to the afore mentioned commands of MEDIBUS (see page 38) the Realtime-Extension uses the following additional commands.

Command Codes


Command	Code
Request Realtime Configurati on	53H
Configure Realtime Transmission	54H
Realtime Configuration changed	56H

Request Realtime Configuration Command

This command is sent in order to request the supported realtime data of the other device.

Configure Realtime Transmission Command

Ahead of this command the "Request Realtime Configuration" Command must be sent to request the supported realtime data. The "Configure Realtime-Transmission" Command is sent to request realtime-data in the order and interval as given in the command's argument. The argument has the following format:



DATA CODE: Two Byte ASCII HEX number specifying the requested realtime-data (refer to Appendix 1 for data-codes).

MULTIPLIER: Two Byte ASCII HEX number defining the multiples of the sample interval of the respective parameter as reported within the response of the "Request Realtime- Configuration"-command before. (A value of 1 causes the device to send each sampled data. If the value is set to 2, every second sample is transmitted, and so on.)

NOTE

Up to 12 r ealtime-data may be configured by MEDIBUS. But due to technical restrictions some devices allow less realtime-data to be configured.


Realtime Configuration Changed Command

This command is sent, if during runtime the configuration has been changed. E.g., the selection of an anaesthesia gas has been changed. The receiving device can now request the current realtime configuration.

MEDIBUS Realtime-Extension Responses

Realtime Configuration Response

This response is sent on reply to the receipt of a Request Realtime-Configuration Command. The response contains the codes of realtime data, that can be sent by the requested device, the related sample interval, minimal and maximal values:



DATA CODE: Two Byte ASCII HEX number which uniquely identifies the realtime data.

INTERVAL: 8 Byte ASCII value specifying the sample interval of the realtime data in ms.

MIN: 5 Byte ASCII minimal decimal value of the realtime data.

NOTE

MIN must be transmitted as a binary 0. Thus all binary data are signless (000hex = 0dec to FFFhex = 4095dec).

MAX MIN must be transmitted as a binary 0. Thus all binary data are signless (000hex = 0dec to FFFhex = 4095dec).

MAXBIN 3 Byte ASCII HEX representation of the binary value which is equal to the maximal-realtime data given in MAX. $MAXBIN = (MAX - MIN)/Resolution$

NOTE

In MIN and MAX the first character may be a minus. MIN and MAX may contain a decimal point at any position except at the first or last character position. Leading zeros must be replaced by spaces. Due to the definitions above the receiver can calculate the value $v(x_{bin})$ of a received binary realtime data x_{bin} by the following equation:

$$v(x_{bin}) = MIN + x_{bin} * (MAX - MIN)/MAXBIN$$


If the data source has an internal binary representation of the realtime data where the binary value x_{MIN} of MIN is not 0 then it can calculate the binary value to transmit x_{bin} from any internal value $x_{internal}$ as follows:

$$X_{bin} = X_{internal} - X_{MIN}$$

MEDIBUS Realtime-Extension Realtime-Data Records

Structure of Realtime-Data Records

A realtime data record has the following format:



022

One Sync byte is transmitted first followed by 0 up to 16 Sync-Commands and up to 12 realtime data.

Sync Byte

The sync byte status bits, if set, hold the following information:

Bit no.	Message
0	Value of 1st Realtime Curve is following
1	Value of 2nd Realtime Curve is following
2	Value of 3rd Realtime Curve is following
3	Value of 4th Realtime Curve is following

A status bit of the Sync Byte is set for each following realtime value. The bit-positions are related to the datastream-number in the same order as specified within the "Configure Realtime-Transmission" command by the data-requesting device before. A missing value is replaced by the value of the next curve. In this case the total number of transmitted data bytes is reduced by 2 or multiple of 2, respectively.

Example (binary):

11010101	10110000 10100011	10001011 10111100
Sync Byte	value of 1st curve	value of 3rd curve
	8F0H	F0BH

2nd value
missing

If more than 4 realtime data have been configured concurrently, their appearance is indicated using the Sync-Commands "Transmitted Datastreams". These commands are transmitted once each time the order or number of the transmitted datastreams 5 to 12 **changes** (e.g. if a data is newly available or missing).

The Sync Byte is sent only heading a data record, never alone.


Sync-Commands

Sync-Commands are used for information which have to be sent very fast (e.g. enable/disable datastreams) or which belong to the realtime data directly following the Sync-Commands.

A Sync-Command is not responded to!

Sync-Commands are embedded within a realtime-data record. Up to 16 Sync-Commands may optionally follow the Sync Byte.

A Sync-Command consists of two Sync-Command-Bytes:



098

The first Sync-Command byte contains the Sync Command Code in bits 0 to 3. The second byte holds the argument in bits 0 to 3 also.

The following Sync-Commands are defined currently

Command Code	Argument	Meaning
1100 0000	1100 0000	End of Sync-Command Sequence
1100 0001	1100 xxxx	Datastream enable/disable (data 1-4) ¹⁾
1100 0010	1100 xxxx	Datastream enable/disable (data 5-8) ¹⁾
1100 0011	1100 xxxx	Datastream enable/disable (data 9-12) ¹⁾
1100 0100	1100 xxxx	Transmitted Datastreams (data 5-8) ^{1) 2)}
1100 0101	1100 xxxx	Transmitted Datastreams (data 9-12) ^{1) 2)}
1100 0110	1100 xxxx	Synchronisation Commands:
	0000	Start of Ventilator Inspiratory Cycle
	0001	Start of Ventilator Expiratory Cycle
1100 1111	1100 0000	Corrupt Data-Record received

- 1) The bits 0 to 3 within the argument are related to the datastreams in the same order as the Sync-Byte's status-bits to the datastream 1 to 4 (see page 17). A set bit (1) enables, a resetted bit (0) disables the attached datastream.
- 2) "Transmitted Datastreams" commands must be sent only if the order of transmitted datastreams has changed. The datastreams will be sent in a fixed order, until a new "Transmitted Datastream" command announces a new order.

Framing of sync-command sequences:

Sync-Commands are always headed by a Sync Byte as described in section 5.1.

Within a data record the end of a sequence of Sync-Commands is indicated by the first Data-Value Byte.

If the Sync-Commands are transmitted without any realtime-data, the Sync Byte (see page 31) must be set to 1101 0000 (no data) and the last Sync-Command Byte must be followed by the Sync-Command "End of Sync- Command Sequence" (11000000 11000000).

Appendix

MEDIBUS-Life-Cycle

Communications Start-up

After starting the device it must send an ICC-command. This command may be repeated each 3 seconds if no character was received or each 10 seconds if characters have been received but no complete response within this time.

On receipt of an ICC-response a "Request Device-Identification"- command may be sent if required. After having received and processed a correct device identification the device is now in an active protocol state. It is ready to answer Data Request Commands or to send any command on its own.


Communications Time-Out

If no character is received within 3 seconds the communication is assumed to be broken. In case of a 3 seconds time-out, communication must be re-initialized by sending an ICC-command. If communication is in an idle state (no command or response pending) the NOP-command may be sent by any device every 2 seconds to keep communication alive and to avoid a 3 seconds time-out.

Terminating Communication

To terminate communication regular the STOP-command must be sent. After receipt of the STOP-response the device may completely stop the communication. The STOP-command must be sent to avoid an error-state and error-communication on the linked device.

MEDIBUS-Life-Cycle



MEDIBUS Realtime-Extension Life-Cycle

Communications Start-Up

For the general startup sequence please refer to MEDIBUS-life Cycle (page 40).

The first step is to request for the supported real-timedata using the "Request Realtime-Data Configuration" Command.

Data-transmission must be configured by the "Configure Realtime-Transmission" Command, with the argument describing the requested realtime data.

All or single realtime-datastreams can be enabled using the "Enable/Disable Datastream" Commands, depending on how many different realtime data have been configured and which shall be enabled.

The enable/disable state of the datastreams can be changed by sending a new "Enable/Disable Datastream" Command with new arguments.

After receipt of an "Enable/Disable Datastream" Command the device starts transmission of the defined realtime data at the defined rate. Transmission is continued until receipt of new "Enable/Disable Datastream" Command or a "Configure Realtime- Transmission" Command.

Since there is no feedback about the correct receipt of a Sync-Command it is recommended to repeat a "Datastream enable/disable"-Sync-Command in an interval of some seconds. This will avoid misinterpretation of datastreams for a long time if a previous Sync-Command has got lost or been changed during transmission.

Transmission is performed in the format described above.

Re-configuration of data-transmission


Re-configuration by receiver:

After realtime-data transmission has been configured and enabled once, using the "Configure Realtime- Transmission" command a new configuration may be sent at any time. The realtime data source will change the configuration after receipt of the configuration command. To keep the moment of change under control the realtime data receiver is recommended to disable all or only the changed datastreams by using the "Enable/Disable Datastream" commands before sending the new configuration and to enable the datastreams after receipt of the response to its configuration command.

Re-configuration by sender:

If the realtime-data source changes its realtime configuration (e.g. due to a change of resolution of any trace) it has to stop all transmissions of realtime-data to avoid misinterpretations of data by the receiving device. As a second step the source has to send the "Realtime Configuration changed" command. Afterwards the receiving device may repeat the realtime-data configuration sequence with a "Request Realtime-Data Configuration" command.

MEDIBUS Realtime-Extension Life-Cycle



ASCII HEX Format

Alarm codes, data codes and checksums are transmitted in ASCII HEX format.

This format uses two bytes to represent a one-byte hexadecimal value in ASCII characters. The HEX numbers A – F are represented in upper case letters only.

Leading zeroes are not suppressed.

Example:


Decimal => Binary => HEX(H) => ASCII HEX Representation

58 => 0011 1010 => 3AH

ASCII "3" (33H)	ASCII "A" (41H)
--------------------	--------------------

MEDIBUS example 1

Communication Start-up



Dräger Device issues the ICC command:

ESC	ICC	CHECKSUM		CR
1BH	51H	36H	43H	0DH

PC responds:

SOH	ICC-Echo	CHECKSUM		CR
01H	51H	35H	32H	0DH

PC issues the device **ID** command:

ESC	Request Device Identity	CHECKSUM		CR
1BH	52H	36H	44H	0DH

Dräger Device responds:

SOH	Command Echo	ID Number	Device Name	Device Revision	MEDIBUS Revision	CHECKSUM		CR
01H	52H	8888	'Draeger Device '	03.00	03.00	38H	43H	0DH

Dräger Device issues the device **ID** command:

ESC	Request Device Identity	CHECKSUM		CR
1BH	52H	36H	44H	0DH

PC responds either with a complete Device Identification:

SOH	ID Echo	ID Number	Device Name	Device Revision	MEDIBUS Revision	CHECKSUM		CR
01H	52H	0161	'My PC'	01.00	03.00	46H	43H	0DH

or with an empty Device Identification:

SOH	ID Echo	CHECKSUM		CR
01H	52H	35H	33H	0DH

During a period of idleness: Dräger Device issues NOP commands every 2 seconds.

ESC	NOP	CHECKSUM		CR
1BH	30H	34H	42H	0DH

MEDIBUS example 2

Data update command and response

To request data from a device, the following command is sent:

ESC	Command CODE	CHECKSUM		CR
1BH	24H	33H	46H	0DH
0	1	2	4	5 Byte

120

If the receiving device has no data available, the response will be:

SOH	Command ECHO	CHECKSUM		CR
01H	24H	32H	35H	0DH
0	1	2	4	5 Byte

121

If the receiving device has »O2 SAT 98 %« and »OXI PULSE 70/min« to send, the response will be:

SOH	Echo	Data Code		O2 SAT ASCII Data 98 %				Data Code		PULSE ASCII Data 70				Checksum		CR
-	-	'E'	'B'	' '	'9'	'8'	' '	'E'	'1'	' '	'7'	'0'	' '	'7'	'A'	-
01H	24H	45H	42H	20H	39H	38H	20H	45H	31H	20H	37H	30H	20H	37H	41H	0DH
0	1	2	4	8	10	14	16	17								

122

MEDIBUS example 3

Configure Data response command

If the receiving device of example 2 is only interested in O₂Sat, it has to send the configuration Data Response Command.

ESC	Command Code	Datatype	Data Code		CHECKSUM		CR
			'E'	'B'	'1'	'0'	
1BH	4AH	24H	45H	42H	31H	30H	0DH

130

The Dräger Device responds:

SOH	Command ECHO	CHECKSUM	CR
01H	4AH	34H 42H	0DH

131

After sending again the request current data command

ESC	Command CODE	CHECKSUM	CR
1BH	24H	33H 46H	0DH

0 1 2 4 5 Byte

132

the receiving device sends:

SOH	Echo	Data Code		O ₂ SAT ASCII Data 98 %				Checksum		CR
-	-	'E'	'B'	' '	'9'	'8'	' '	'7'	'A'	-
01H	24H	45H	42H	20H	39H	38H	20H	37H	41H	0DH

0 1 2 4 8 11 Byte

133

MEDIBUS example 4

Request current date and time

After Start-up a PC requests current date and time from Dräger Device:

ESC	Request current Date and Time	CHECKSUM		CR
1BH	28H	34H	33H	0D

140

The Dräger Device responds as follows:

SOH	Echo	Time										Date						Checksum		CR	
-	-	'1'	'1'	'.'	'1'	'1'	'.'	'1'	'9'	'1'	'1'	'.'	'N'	'O'	'V'	'.'	'9'	'3'	'E'	'6'	
01H	28H	31H	31H	3AH	31H	31H	3AH	31H	39H	31H	31H	2DH	4EH	4FH	56H	2DH	39H	33H	45H	36H	0DH

141

After a while the user changes the time and date setting on the Dräger Device, which therefore sends the "Time changed" command:

ESC	Time changed	CHECKSUM		CR
1BH	49H	36H	34H	0DH

142

The PC echos

SOH	Echo	CHECKSUM		CR
01H	49H	34H	41H	0DH

143

After reception of the "Time changed" command the PC is informed about a new time and date setting and may now, before, while or after transmitting the command-echo send again the "Request current date and time" command to get the new setting.

MEDIBUS example 5

Trend Data Transmission

A PC wants to request trend data from a Dräger device.
First, it requests the actual trend data status:

ESC	Request Trend Data Status	CHECKSUM	CR
1BH	6CH	38H 37H	0D

The Dräger device supports the trend parameters "Derived Respiratory Rate" (code page 1, data code D9H) and "Oxygen Saturation" (code page 1, data code EBH).

In this example 2345 (= 929_{hex}) Respiratory Rate samples are available, in which the eldest sample was generated on 16.08.2003 at 12.34:56 (time stamp = 2F10631C_{hex} with YYYYYY = 23, MMMM = 8, DDDD = 16, hhhh = 12, mmmmm = 24, ssss = 56).

There are 912 (= 390_{hex}) Oxygen Saturation samples starting on 16.08.2003 at 13.05:10 (time stamp = 2F1068A5_{hex}).

SOH	Echo	Parameter 2: Derived Respiratory Rate																	
		Code Page	Data Code		Count				Begin										
		'2'	'4'	'D'	'9'	''	''	''	'9'	'2'	'9'	'2'	'F'	'1'	'0'	'6'	'3'	'1'	'C'
01H	6CH	32H	34H	44H	39H	20H	20H	20H	39H	32H	39H	32H	46H	31H	30H	36H	33H	31H	43H

Code Page	Data Code	Count														Begin				Checksum	CR
'2'	'4'	'E'	'B'	''	''	'3'	'9'	'0'	'2'	'F'	'1'	'0'	'6'	'8'	'A'	'5'	'B'	'0'			
32H	34H	45H	42H	20H	20H	20H	33H	39H	30H	32H	46H	31H	30H	36H	38H	41H	35H	42H	30H	0DH	

The PC requests all Respiratory Rate samples, beginning with the eldest available sample in sequences of 100 samples. Therefore it starts with a request for 100 (= 64_{hex}) trend data, using the time stamp of the eldest sample:

ESC	Request Trend Data	Code Page	Data Code	Count	Begin								Checksum	CR				
		'2'	'4'	'D'	'9'	'6'	'4'	'2'	'F'	'1'	'0'	'6'	'3'	'1'	'C'	'9'	'B'	
1BH	6DH	32H	34H	44H	39H	36H	44H	32H	46H	31H	30H	36H	33H	31H	43H	39H	42H	0DH

The Dräger device sends 100 Respiratory Rate samples, beginning with the eldest requested trend sample:

SOH	Echo	Code Page	Data Code	Count	Value 1				Time Stamp 1										
		'2'	'4'	'D'	'9'	'6'	'4'	'1'	'4'	''	''	'2'	'F'	'1'	'0'	'6'	'3'	'1'	'C'
01H	6DH	32H	34H	44H	39H	36H	44H	31H	34H	20H	20H	32H	46H	31H	30H	36H	33H	31H	43H
.....value 2 to value 99.....																			
		Value 100				Time Stamp 100				Checksum	CR								
		'1'	'6'	''	''	'2'	'F'	'1'	'0'	'6'	'5'	'5'	'C'	'?	'?	'0DH			
31H	36H	20H	20H	32H	46H	31H	30H	36H	35H	35H	43H								

In order to get the next 100 trend samples, the PC sends another trend data request. To achieve a correct successive row of trend samples, the PC uses the time stamp of the last sample in the preceding trend data sequence, incremented by one.

ESC	Request Trend Data	Code Page	Data Code	Count	Begin								Checksum	CR				
		'2'	'4'	'D'	'9'	'6'	'4'	'2'	'F'	'1'	'0'	'6'	'5'	'5'	'D'	'A'	'2'	
1BH	6DH	32H	34H	44H	39H	36H	44H	32H	46H	31H	30H	36H	35H	35H	44H	41H	32H	0DH

The Dräger device will respond with the next successive 100 trend samples.

This is repeated until all available trend samples have been transmitted to the PC.

MEDIBUS example 6

Realtime-Extension

After general start-up procedure a PC requests a Dräger Device for its Realtime Configuration:

ESC	Request Realtime Configuration	CHECKSUM		CR
1BH	53H	36H	45H	0D

210

The Dräger Device offers Airway Pressure (Code 00H) in a range from –10 to 100 mbar with a resolution of 1/8 mbar (MAXBIN = (MAX – MIN)/Resolution = 880_{dec} = 370_{hex}). And CO₂ (Code 06H) in a range from –20 to 100 mmHg, where the MAX-value 100 mmHg is transmitted as MAXBIN = 1120_{dec} = 460_{hex}. Both values will be transmitted at a minimal interval of 16 ms.

SOH	Echo
01H	53H

Code		Interval									Minimal Value				Maximal Value				MaxBin				
'0'	'0'	' '	' '	' '	'1'	'6'	'0'	'0'	'0'	'0'	'-	' '	' '	'1'	'0'	' '	' '	'1'	'0'	'0'	'3'	'7'	'0'
30H	30H	20H	20H	20H	31H	36H	30H	30H	30H	30H	2DH	20H	20H	31H	30H	20H	20H	31H	30H	30H	33H	37H	30H

Code		Interval									Minimal Value				Maximal Value				MaxBin			
'0'	'6'	' '	' '	' '	'1'	'6'	'0'	'0'	'0'	'-	' '	' '	'2'	'0'	' '	' '	'1'	'0'	'0'	'4'	'6'	'0'
30H	36H	20H	20H	20H	31H	36H	30H	30H	30H	2DH	20H	20H	32H	30H	20H	20H	31H	30H	30H	34H	36H	30H

Checksum		CR
'3'	'B'	
33H	45H	0DH

211

PC now requests each second sample of Airway Pressure and each third sample of CO₂, respectively:

ESC	Configure Realtime Transmission Command Code	Code		Multipl		Code		Multipl		Checksum		CR
		'0'	'0'	'0'	'2'	'0'	'6'	'0'	'3'	'F'	'A'	
1BH	54H	30H	30H	30H	32H	30H	36H	30H	33H	46H	41H	0DH

212

Dräger Device sends a command acception response without transmitting any realtime data yet:

SOH	Echo	Checksum	CR
01H	54H	35H	0DH

213

PC enables both realtime traces:.

Sync Byte	Sync Command enable Traces 1 + 2		End of Sync Command Sequence	
D0H	C1H	C3H	C0H	C0H

214

Dräger Device starts periodical transmission of Realtime Data:

time $t = 0$:

One realtime value for each trace is sent:

Sync Byte	Realtime Value Airway Pressure		Realtime Value CO ₂	
D3H	91H	81H	8DH	83H

Sync Byte =

D3H = 1101 0011_{bin} indicates realtime values 0 and 1 are following because bits 0 (LSB) and 1 are set.

Airway Pressure =

91H 81H = 1001 0001_{bin} 1000 0001_{bin}

=> Regarding the format of Realtime-Value-Bytes (see page 32) the twelve bit realtime value is 0000 0101 0001_{bin} = 81_{dec}

=> the according Airway Pressure value $v(x_{bin})$ in mbar is (regard page 35):

$$v(x_{bin}) = -10 \text{ mbar} + 81 * (100 \text{ mbar} - (-10 \text{ mbar})) / 880 = 0.125 \text{ mbar}$$

CO₂ =

8DH 83H = 1000 1101_{bin} 1000 0011_{bin}

=> twelve bit binary value 0000 1100 1101_{bin} = 205_{dec}

=> CO₂ value $v(x_{bin})$ in mmHg:

$$v(x_{bin}) = -20 \text{ mmHg} + 205 * (100 \text{ mmHg} - (-20 \text{ mmHg})) / 1120 = 1.96 \text{ mmHg}$$

time $t = 0 + 1 * 16 \text{ ms}$:

No Realtime values are sent

time $t = 0 + 2 * 16$ ms:

Only a value for Airway Pressure is sent:

Sync Byte	Realtime Value Airway Pressure	
D1H	94H	81H

216

Sync Byte D1H = 1101 0001_{bin} has only bit 0 set, indicating only the first configured trace Airway Pressure is following the Sync Byte. The second trace CO₂ is missing.

time $t = 0 + 3 * 16$ ms:

Only a value for CO₂ is sent:

Sync Byte	Realtime Value CO ₂	
D2H	8EH	83H

220

Sync Byte D2H = 1101 0010_{bin} has only bit 1 set, indicating only the second configured trace CO₂ is following the Sync Byte. The first trace Airway Pressure is missing.

At each start of ventilator inspiratory cycle Dräger Device may inform PC about that event so PC may synchronize its Airway Pressure displaying. In this case Dräger Device will send the "Start of ventilator inspiratory cycle" Sync-Command with the respective realtime value for Airway Pressure:

Sync Byte	Start of Insp. Cycle		Realtime Value Airway Pressure		Realtime Value CO ₂	
D3H	C6H	C0H	90H	81H	8AH	83H

221

Frequently Asked Questions

The following are frequently asked questions showing some typical problems which occur when, for instance, a PC is being programmed to communicate with a Dräger medical device. The respective answers should help the understanding of and provide advice for the programming process.

NOTE

Within the MEDIBUS protocol definition Dräger medical devices may show varying behavior. For this reason, the following general answers may not apply in each and every case for each and every device.

Do all Dräger medical devices send the ICC command on their own?

Yes. As long as an ICC response and an ICC command are not received, Dräger medical devices send the ICC command approx. every 3 seconds.

Does my PC have to send an ICC command itself, or will it suffice to respond to an ICC command sent by a Dräger medical device?

Responding to the ICC command sent by the Dräger medical device will suffice.

Does my PC have to respond to the ICC command from the Dräger medical device?

Yes. Every command must be responded to.

What kind of identification should my PC send?

For most Dräger medical devices, it is enough to send an empty response. Some devices (e.g. Cicero B and Evita 4) do, however, expect a complete response with specific content. For this reason, it is recommended that a complete response containing the following information is sent:

ID Number:	0161
Name of Device:	any
Device Version:	any

MEDIBUS Version: the version number for the MEDIBUS protocol definition according to which the PC program was developed

Compare with MEDIBUS Example 1 in this protocol definition.

Does my PC also have to send a device identification request?

Only if your PC program requires this information, otherwise: No

Why isn't the data configuration command processed; in other words: Why does my PC still receive all data when subsequent data request commands are given?

The data configuration command is often sent too early. All Dräger medical devices send a device identification request after the ICC sequence. Your PC needs to respond to this request first (see above) and then allow the Dräger medical device approx. 200 ms to process this response. Only after this will the Dräger medical device be ready to process the data configuration command. All previous data configuration commands shall be ignored by the Dräger medical device.

My PC sets up the connection to the Dräger medical device and communicates for some time with it. Why does the Dräger medical device send another ICC command after some time?

In many cases, the command sent by the Dräger medical device is answered either incorrectly or not at all. A command incorrectly responded to will lead to the command being repeated in a short span of time. A command that is not responded to will be repeated after approx. 10 seconds. After a maximum of 3 repetitions, the Dräger medical device assumes that the connection is faulty and tries to re-establish communication using an ICC command. Thus, if communication only stands for a

few seconds, then your PC is probably sending an incorrect response (e.g. a response with an incorrect format or check sum). If communication is stable for approx. 30 to 40 seconds before the Dräger medical device sends another ICC command, then your PC probably isn't responding to the command from the Dräger medical device at all. (Embedded commands often go unnoticed. In many cases, for example, the device identification request from the Dräger medical device is embedded in the ICC response.)

Will my PC also be asked for data?

If your PC sends a device identification response with the information given above, then: No.

Which commands shall be sent to my PC from the Dräger medical device?

There is, for every device, a table of "Transmitted Commands" in the instructions for use "MEDIBUS for Dräger Anesthesia Devices", "MEDIBUS for Dräger Pediatric Devices" and "MEDIBUS for Dräger Intensive Care Devices".

NOTE

A flexible PC program design that allows every command received from the Dräger medical device to be responded to at any given time has often proven very valuable. In such cases, it should make no difference whether or not the command is known and, as expected, comes at a specific time, or if the command comes at an unexpected time or is even completely unknown. Communicating with the Dräger medical device should pose no problem if every command is responded to with at least an empty response (for exceptions see "device identification response", page 14). Commands which are not responded to, however, always lead to an interruption in communication.

When does my PC has to send or respond to an NOP command?

An NOP command should be sent when other commands (e.g. data request commands) should not be sent to the other device for more

than 2 seconds, when communication should continue to stand and, for this reason, an exceeding of the 3-second timeout has to be prevented.

The Dräger medical device shall send an NOP command if it has received neither command nor response from the PC for approx. 2 seconds, and if it itself did not send a previous command and is, therefore, not waiting for a response. Thus, it will suffice if the PC sends an NOP response to the Dräger medical device within approx. 1 second. It is, however, better and provides better protection against communication errors caused by exceeding time-outs if the PC sends an NOP command itself in such cases.

How often can my PC send a command?

As soon as the response to the preceding command has been received, the next command can be sent. If the PC sends a command before the Dräger medical device could respond in full to the preceding command, the Dräger medical device shall ignore the previous command and begin processing and responding to the new command.

How long does it take the response to come from the Dräger medical device?

This depends on the Dräger medical device, the command, the length of the response and the Baud rate. For this reason, this question cannot be answered completely. Usually, far less than 500 ms pass between the completed reception of a command and the transfer of the first character of the response.

What happens after a 3-second time-out?

If the Dräger medical device detects a 3-second time-out (no signal from the PC for 3 seconds or more), communication shall be re-initialized. Insofar as the Dräger medical device has a corresponding possibility, a communication error will be optically and acoustically displayed until the re-initialization has been successfully completed.

What happens after a 10-second time-out?

When a 10-second time-out occurs (no complete response to a command takes place within 10 seconds), the last command sent by the Dräger medical device shall be repeated; up to three times. After that, communication is re-initialized and the communication error, as was the case for the 3-second time-out, displayed.

Does my PC have to monitor 3-second and 10-second time-outs? Or what happens if these time-outs cannot be monitored or cannot be monitored accurately by my PC?

The 3-second time-out allows the PC to detect whether or not communication with the Dräger medical device stands or has been interrupted. Communication can be re-initialized with an ICC command if there was an interruption and the user would see the communication error displayed. The 10-second time-out allows the PC to monitor whether or not the last command sent to the Dräger medical device was completely and correctly received. In the case of an error, the last command can be repeated or a new command can be sent. If the PC does not have these kinds of requirements, the timeouts on the PC either do not need to be monitored or only monitored less strictly. (Dräger medical devices always monitor these time-outs.)

Is all information and are all configurations sent to the Dräger medical device before the re-initialization lost during a communication re-initialization (ICC), e.g. after a 3-second or a 10-second time-out, and does everything have to be re-sent for this reason?

Yes. All information in the Dräger medical device sent previously will be deleted with every communication (re-)initialization.

Where can commands be embedded in responses?

Commands can be embedded in responses at any position.

Where can real-time data records be embedded in responses and commands?

Real-time data records can be embedded in responses and commands at any position.

Will the ESC for commands / the SOH for responses be calculated into the check sum?

Yes. See the MEDIBUS examples in this instruction for use.

What influence does a command which is embedded in a response have on the check sum calculation?

None. Check sums from commands and responses are always calculated separately. Thus, a command embedded in a response does not influence the check sum of the response or vice versa.

Do real-time data records embedded in commands or responses influence the check sum?

Not at all. Real-time data records are ignored by the check sum calculation.

What does the Stop command do?

Sending a Stop command to the Dräger medical device prevents the Dräger medical device from constantly displaying communication errors after the communication has been ended.

Logbook of Changes


Version	Changes	backward compatible
2.00	First released version	- - -
3.00	Protocol extended by <ul style="list-style-type: none"> - Device Settings - Text Messages - Realtime transmission 	yes
4.00	Protocol extended by second codepage for current measured Data, Alarm Limits and Alarms	yes
4.01	New chapter "Frequently Asked Questions".	yes
	No changes of protocol	
4.02	Correction of chapter	no
	"Time & Date Update Response"	
4.03	New Realtime Sync-Command "Start of Ventilator Expiratory Cycle"	yes
5.0	Protocol extended by Trend Data	yes
6.0	Protocol extended by 3rd alarm code page	yes
	Extended device protocol for syringe pumps	

Table of used ASCII-Codes

Code	Explanation	Hexadecimal code
SOH	Start of Header	01H
	Start of Respondings	
ETX	End of Text	03H
CR	Carriage Return	0DH
	End of Commands or Respondings	
DC1	Device Control 1	11H
	Transmitting Device suspends Transmission immediately	
DC3	Device Control 3	13H
	Transmitting Device resume Transmission	
NAK	Negative Acknowledge	15H
	Character for corrupt Commands	
CAN	Any Transmission in Progress will be immediately aborted.	18H
ESC	Escape	1BH
	Start of Commands	
Space	Blank Key	20H
	Surplus Character Positions (Data and Settings) must be filled up with this Character.	
	Character to fill up leading Zeroes.	
20H – 7FH	Range used for the "slow" MEDIBUS Protocol	
80H – FFH	Range used for the MEDIBUS Realtime-Extension	



Directive 93/42/EEC
concerning Medical Device

 Manufacturer

 **Dräger Medical AG & Co. KG**
Moislinger Allee 53 – 55
D-23542 Lübeck
Germany
 +49 451 8 82-0

FAX +49 451 8 82-20 80
 <http://www.draeger.com>

90 28 258 – TD 6494.380 en
© Dräger Medical AG & Co. KG
10th edition – December 2007
Dräger Medical reserves the right to make
modifications to the equipment without prior
notice.