


# Alaris® Pump module, model 8100 pressure sensors tip sheet

This tip sheet provides information on the large volume pump (LVP) membrane and pressure sensors.



The two pressure sensors are covered by a clear, rectangular membrane shown in the picture at left.

**CAUTION:** To avoid pressure sensor damage, do not apply pressure to the center of the pressure sensors.

If you have questions about the Alaris® System, please visit our infusion documentation website at [carefusion.com/InfusionDocumentation](http://carefusion.com/InfusionDocumentation) or contact us at [customerfeedback@carefusion.com](mailto:customerfeedback@carefusion.com).

© 2015 CareFusion Corporation or one of its affiliates. All rights reserved. Alaris, CareFusion and the CareFusion logo are trademarks or registered trademarks of CareFusion Corporation or one of its affiliates. IF5174 (0715/250)

CareFusion  
San Diego, CA

[carefusion.com](http://carefusion.com)

