











**Attachment E: Ordering Supplies - Old vs. New**





**IV Sets and Accessories Catalog - Dedicated**



**Updated: April 19, 2016**










<b><u>Current Storeroom #</u></b>	<b><u>Current Product</u></b>	<b><u>New Storeroom #</u></b>	<b><u>UVA Description</u></b>	<b><u>New Product</u></b>	<b><u>Areas to be stocked</u></b>
<p><b><u>97240</u></b> <b><u>93271</u></b></p>	<p><b><u>72023E</u></b></p>  <p>zavlev seft-alkeson S, pvlav kzerb Jber nolavfml kcol teul slem esasp-5 mont "E bns" EB kroqps mo A2S "001" rrtqpsl fee slemkroqpsA Jim SS semulov priming slemkroqpsA</p>	<p><b><u>97250</u></b></p>	<p><b><u>Tubing IV Pump Primary</u></b></p>	<p><b><u>2426-0007 (Primary)</u></b></p>  <p>Alaris Pump Module set Check Valve Roller, Slide Clamp(s) <b>3 SmartSite needle-free valve(s)</b> 6", 74", 90" (1 above pumping segment and 2 below) from 2-piece Male Luer Lock. Not Made with DEHP. 20 Drop L: 126 in L: 320 cm PV: 26 mL Fluid Path Sterile</p>	<p><b><u>Everywhere</u></b></p>
<p><b><u>97241</u></b> <b><u>93274</u></b></p>	<p><b><u>72980E</u></b></p>  <p>Non-vented blood set, 180 micron filter, 1 needle-free valve 9" from 2-piece male luer lock Approximate set length: 105" 267 cm. approx. Approximate priming volume: 31 mL</p>	<p><b><u>97251</u></b></p>	<p><b><u>Tubing IV Pump Blood Infusion</u></b></p>	<p><b><u>2477-0007</u></b></p>  <p>Alaris Pump Module set Non-Vented <b>Blood Set</b> 180 Micron Filter Roller, 2 Slide Clamp(s) 1 SmartSite needle-free valve(s) 65" (below pumping segment) from 2-piece Male Luer Lock. Not Made with DEHP. 15 Drop L: 121 in L: 308 cm PV: 33 mL Fluid Path Sterile</p>	<p><b><u>Inpatient Anesthesia</u></b></p> <p><b><u>*Do NOT stock in NICU</u></b></p>
<p><b><u>97243</u></b> <b><u>99415</u></b></p>	<p><b><u>72953</u></b></p>  <p>Infusion set, low sorbing tubing segment, 2 piece male luer lock Approximate set length: 110" 279 cm. approx. Approximate priming volume: 22 mL</p>	<p><b><u>97257</u></b></p>	<p><b><u>Tubing Low Sorbing Short Set</u></b></p>	<p><b><u>11426964 (Low Sorb Short Chemo Set for Oncology Areas ONLY)</u></b></p>  <p>Alaris Pump Module set <b>Low Sorbing</b> Roller Clamp(s) 2-piece Male Luer Lock. Not Made with DEHP. 20 Drop L: 47 in L: 120 cm PV: 15 mL Fluid Path Sterile</p>	<p><b><u>Oncology Areas</u></b></p> <p><b><u>ECCC</u></b></p> <p><b><u>Battle</u></b></p> <p><b><u>7 Acute</u></b></p> <p><b><u>8 West</u></b></p> <p><b><u>Stem Cell</u></b></p> <p><b><u>Hope</u></b></p>

<p><b><u>97242</u></b></p>	<p style="text-align: center;"><b><u>72923</u></b></p>  <p>Infusion set, low sorbing tubing segment, 0.2 micron (low-protein binding) filter, 2-piece male luer lock  Approximate set length: 104" 264 cm, approx.  Approximate priming volume: 24 mL</p>	<p><b><u>97258</u></b></p>	<p><b><u>Tubing Low Sorb EXT 0.2 Filter</u></b></p>	<p style="text-align: center;"><b><u>20350E</u></b>  <b>(Add On Low Sorb Extension Set w/ 0.2 Micron Filter for Oncology Areas ONLY)</b></p>  <p><b>Extension set Low Sorbing 0.2 Micron Filter</b> (1.72 sq. in Low Protein Binding) Pinch Clamp(s) 1 SmartSite needle-free valve(s) 6" from Spin Male Luer Lock. Not Made with DEHP. L: 17 in L: 43 cm PV: 5 mL Fluid Path Sterile</p>	<p><b><u>Oncology Areas</u></b>  <b><u>ECCC</u></b>  <b><u>Battle</u></b>  <b><u>7 Acute</u></b>  <b><u>8 West</u></b>  <b><u>Stem Cell</u></b>    <b><u>Hope</u></b></p>
<p><b><u>97244</u></b></p>	<p style="text-align: center;"><b>Currently using Secondary Set</b></p> <p style="text-align: center;"><b>Codan</b></p>	<p><b><u>97261</u></b></p>	<p><b><u>Tubing IV Pump Secondary</u></b></p> <p>(Description update only)</p>	<p style="text-align: center;"><b><u>2204-0007</u></b>  <b>(Short Chemo Set for Oncology Areas ONLY)</b></p>  <p>Alaris Pump Module set Roller Clamp(s) 2-piece Male Luer Lock. Not Made with DEHP. 20 Drop L: 47 in L: 119 cm PV: 14 mL Fluid Path Sterile</p>	<p style="text-align: center;"><b><u>Everywhere</u></b></p>
<p><b><u>99416</u></b> <b><u>97242</u></b></p>	<p style="text-align: center;"><b><u>72923</u></b></p>  <p>Infusion set, low sorbing tubing segment, 0.2 micron (low-protein binding) filter, 2-piece male luer lock  Approximate set length: 104" 264 cm, approx.  Approximate priming volume: 24 mL</p>	<p><b><u>97259</u></b></p>	<p><b><u>Tubing Low Sorbing Primary Set</u></b></p>	<p style="text-align: center;"><b><u>24600-0007</u></b>  <b>(Low Sorb 3-Port Set for Oncology Areas ONLY)</b></p> <p style="text-align: center;"><b>No Image Available</b></p> <p>Alaris Pump Module set <b>Low Sorbing</b> Check Valve Roller Clamp(s) <b>3 SmartSite needle-free valve(s)</b> 6", 66", 88" (1 above pumping segment and 2 below) from 2-piece Male Luer Lock. Not Made with DEHP. 20 Drop L: 120 in L: 304 cm PV: 25 mL Fluid Path Sterile</p>	<p><b><u>Oncology Areas</u></b>  <b><u>ECCC</u></b>  <b><u>Battle</u></b>  <b><u>7 Acute</u></b>  <b><u>8 West</u></b>  <b><u>Stem Cell</u></b>    <b><u>Hope</u></b></p>

<p><b><u>93272</u></b></p> <p><b><u>Delete from Stock</u></b></p>	<p><b><u>72103E</u></b></p>  <p>Infusion set, 150 mL, burette (needle-free valve), 1 needle-free valve 9" from 2-piece male luer lock Approximate set length: 120" 305 cm, approx. Approximate priming volume: 22 mL</p>			<p><b>Removed from Stock</b></p>	
<p><b>Replacing</b></p> <p><b><u>93293</u></b></p> <p><b><u>93295</u></b></p>	<p><b>N/A</b></p>	<p><b><u>97252</u></b></p>	<p><b><u>Tubing IV Pump Narcotic No Port</u></b></p>	<p><b><u>2200-0500</u></b> <b>(No Port Set for Use w/ Bag Based Narcotics)</b></p>  <p>Alaris Pump Module set Roller Clamp(s) 2-piece Male Luer Lock. 20 Drop L: 109 in L: 277 cm PV: 23 mL Fluid Path Sterile</p>	<p><b><u>Inpatient and Procedural areas</u></b></p>
	<p><b><u>N/A</u></b></p>	<p><b><u>97253</u></b></p>	<p><b><u>Tubing IV Pump Syringe PSD</u></b></p>	<p><b><u>10014914</u></b> <b>(Pressure Sensing Disc for NICU/Anesthesia ONLY)</b></p>  <p>Alaris Syringe Module set Microbore <b>Pressure Sensing Disc</b> Slide Clamp(s) Fixed Male Luer Lock. Not Made with DEHP. L: 60 in L: 152 cm PV: 0.7 mL Fluid Path Sterile</p>	<p><b><u>NICU and Anesthesia only</u></b></p>
	<p><b>N/A</b></p>	<p><b><u>97254</u></b></p>	<p><b><u>Tubing IV Pump Syringe 0.2 PSD</u></b></p>	<p><b><u>10014916</u></b> <b>(Pressure Sensing Disc w/ 0.2 Micron Filter for NICU/Anesthesia ONLY)</b></p>  <p>Alaris Syringe Module set Microbore <b>Pressure Sensing Disc 0.2 Micron Small Volume Filter</b> (0.25 sq. in Low Protein Binding) Slide Clamp(s) Fixed Male Luer Lock. Not Made with DEHP. L: 60 in L: 152 cm PV: 1 mL Fluid Path Sterile</p>	<p><b><u>NICU and Anesthesia only</u></b></p>

N/A	N/A	<u>97260</u>	<u>Tubing Low Sorbing Syringe Set</u>	<p><b><u>C20014</u></b>  <b>(For 8W/Stem Cell ONLY)</b></p>  <p>Extension set Low Sorbing 2-piece Male Luer Lock. Not Made with DEHP. L: 30 in L: 76 cm PV: 4 mL Fluid Path Sterile</p>	<u>8W Stem cell</u>
N/A	N/A	<u>97255</u>	<u>Tubing IV Pump Syringe Adaptor</u>	<p><b><u>10010483</u></b>  <b>(Syringe Adapter Set for Syringe Infusion via LVP – for Anesthesia ONLY)</b></p>  <p>Alaris Pump Module set <b>Syringe Adapter</b> Smallbore Slide Clamp(s) 2-piece Male Luer Lock. Not Made with DEHP. L: 82 in L: 208 cm PV: 4 mL Fluid Path Sterile</p>	<u>Anesthesia</u>

<b>Alaris PCA Module</b>					
<b><u>Current Storeroom #</u></b>	<b><u>Current Product</u></b>	<b><u>New Storeroom #</u></b>	<b><u>UVA Description</u></b>	<b><u>New Product</u></b>	<b><u>Areas to be Stocked</u></b>
<b><u>93293</u></b>	<p style="text-align: center;"><b><u>21-7047-24</u></b></p>  <p>CADD® Extension Set with male luer, clamp, and integral anti-siphon valve with male luer 60"</p>	<b><u>97256</u></b>	<b><u>Tubing PCA IV Pump</u></b>	<p style="text-align: center;"><b><u>30873</u></b></p>  <p>Alaris PCA Module set Microbore Check Valve Anti-Siphon Valve 2 Pinch Clamp(s) Y-Connector 2 Female Luers 7" from Fixed Male Luer Lock. Not Made with DEHP. L: 90 in L: 229 cm PV: 3 mL Fluid Path Sterile</p>	<b><u>Everywhere</u></b>
<b><u>93295</u></b>	<p style="text-align: center;"><b><u>21-7092-24</u></b></p>  <p>CADD® Extension Set with male luer, integral anti-siphon valve with male luer, y-extension with one-way female checkvalve, and clamps 68"</p>	<b><u>97256</u></b>	<b><u>Tubing PCA IV Pump</u></b>	<p style="text-align: center;"><b><u>30873</u></b></p>  <p>Alaris PCA Module set Microbore Check Valve Anti-Siphon Valve 2 Pinch Clamp(s) Y-Connector 2 Female Luers 7" from Fixed Male Luer Lock. Not Made with DEHP. L: 90 in L: 229 cm PV: 3 mL Fluid Path Sterile</p>	<b><u>Everywhere</u></b>
<b><u>Not currently carried at UVA</u></b>	<p style="text-align: center;"><b><u>30843E</u></b></p>  <p>PCA extension set, smallbore tubing segment, check valve, 1 needle-free valve 7" from 2-piece male luer lock</p> <p>Approximate set length: 13" 33 cm. approx.</p> <p>Approximate priming volume: 1 mL</p> <div style="text-align: right;">  </div>			<p style="text-align: center;"><b><u>30873</u></b></p>  <p>Alaris PCA Module set Microbore Check Valve Anti-Siphon Valve 2 Pinch Clamp(s) Y-Connector 2 Female Luers 7" from Fixed Male Luer Lock. Not Made with DEHP. L: 90 in L: 229 cm PV: 3 mL Fluid Path Sterile</p>	