


EPIC DOWNTIME

Process for EPIC downtime and initiating a patient on Glucommander:


1. Desktop ICON (Blue “G”) downtime login



2. Ensure the appropriate unit is displayed in the status bar, then select “Add Patient”



3. Select “IV” button



Patient Name	MRN	Date of Birth	Facility Unit	Room
A, TRAINING	M990001017...	06/09/1969	UVHE: 4 NORTH	101
AA, TRAINING	M990001277...	06/09/1969	UVHE: 4 NORTH	127